

Hi,

You will find a manufacturer that is able to make translucent patches that are applied to the human skin.


- The dimensions are roughly 5.2cm by 13cm. Each one patch has 10 circular wells which measure about 0.81cm.
 - These measurements do not have to be exact.
 - Ideally, each of these wells is indented slightly.
 - If you cannot indent the wells then they will have a cotton type circle.
- Each one set consists of two patches
 - That means if I ask for 10 sets it will be 20 patches total.
 - Patch #1 is not the same as Patch #2
- The tap and adhesive used is equivalent to that in non-woven, occlusive, hypoallergenic surgical tape
 - The adhesive is acrylate-based. There is to be no natural rubber, latex, balsams, or rosins in the adhesive or tape.
- Each ingredient is incorporated/ixed by an inactive ingredient which is petrolatum gel (aka "Vaseline")
 - Other inactive ingredients may include hydroxypropyl cellulose, povidone, povidone with butylhydroxytoluene, povidone with sodium bicarbonate and sodium carbonate, and hydroxypropylcellulose with beta-cyclodextrin.
- Please note this is not a request for any illegal narcotics, opioids or controlled and regulated substances.
- You MUST have previously performed a similar job.

What you must do in order to complete the project:

- Deliver and ship 10 sample patch sets without any ingredients for inspection AND then deliver and ship 100 patch sets (= 200 total patches) that do have ingredients.
- You will then supply the name of the manufacturer and the contact person's name and email. If we continue to purchase patches from them you will continue to get a commission.
- You will NOT be paid for providing websites, email addresses, phone numbers of places. A simple Google search can do that.
- You will NOT be paid for 'trying hard' and 'contacting lots of places.' You will not be paid for effort. You WILL be paid for finding a company that can supply a sample patch and patch with ingredients.

Images

- Sketch of the patch



-
- An example of the indented wells:
 - Link:
<https://www.smartpractice.com/shop/wa/style?id=DMAL7000&m=SPA&cid=5082>



- Photo of the indentation with a circular material



- Picture showing indentations



- - *Ideally* the patch would be translucent



The ingredients

Well No.	Compound	CAS	Composition (% w/w)
1	Nickel sulfate	10101-97-0	2.5% pet
2	Cobalt chloride	7791-13-1	1.0% pet
3	Fragrance Mix 1		
	Cinnamyl alcohol	104-54-1	1.0% pet
	Cinnamal	104-55-2	1.0% pet
	Hydroxycitronellal	107-75-5	1.0% pet
	Amyl cinnamal	122-40-7	1.0% pet
	Geraniol	106-24-1	1.0% pet
	Eugenol	97-53-0	1.0% pet
	Isoeugenol	97-54-1	1.0% pet
	Oakmoss absolute	90028-68-5	1.0% pet
	Sorbitan sesquioleate (emulsifier)	8007-43-0	5% pet
4	Balsam of Peru	8007-00-9	25% pet
	Sorbitan sesquioleate (emulsifier)	8007-43-0	5% pet
5	Paraben mix		
	Methylparaben	99-76-3	3% pet
	Ethylparaben	120-47-8	3% pet
	Propylparaben	94-13-3	3% pet
	Butylparaben	94-26-8	3% pet
6	Gluten	8002-80-0	4.5% pet
7	Fragrance Mix 2		
	Hexyl cinnamic aldehyde	101-86-0	5.0% pet

	Coumarin	91-64-5	2.5% pet
	Farnesol	4602-84-0	2.5% pet
	Lylal	31906-04-4	2.5% pet
	Citral	5392-40-5	1.0% pet
	Citronellol	106-22-9	0.5% pet
8	Iodopropynyl butylcarbamate	55406-53-6	0.2% pet
9	Colophonium	8050-09-7	20% pet
10	Quaternium-15	51229-78-8	2.0% pet
11	Corticosteroid mix		
	Hydrocortisone-17-butyrate	13609-67-1	1.0% pet
	Tixocortol-21-pivalate	55560-96-8	1.0% pet
	Budesonide	51333-22-3	0.1% pet
12	PPD	106-50-3	1.0% pet
13	Carba mix		
	1,3-diphenylguanidine	102-06-7	1.0% pet
	Zinc diethyldithiocarbamate	14324-55-1	1.0% pet
	Zinc dibutyldithiocarbamate	136-23-2	1.0% pet
14	Melamine formaldehyde		
	Melamine formaldehyde	9003-08-1	7.0% pet
	Sorbitan sesquioleate (emulsifier)	8007-43-0	1.0% pet
15	Benzophenone-3	131-57-7	10% pet
16	Methyldibromoglutaronitrile	35691-65-7	0.5% pet
17	Lanolin	8027-33-6	30% pet
18	Neomycin sulfate	1405-10-3	20% pet

19	Bacitracin	1405-87-4	20% pet
20	Petrolatum	8009-03-8	100% pet
	-		